

*Thank you and congratulations for purchasing this examination lamp which has been subject of all our care during manufacture and quality control. This product will give you full satisfaction for many years. Before using it, you should read this notice carefully.*

## Legal Compliance:

CE marking according to Directive 93/42 CEE and EN 60598-2-1/A11, EN 60598-2-1, EN 60598-2-4, EN 60598-2-6, EN 60598-2-25, EN 60601-1, EN 60601-1-2, EN 60601-2-41, EN 50581-2-41 (ROHS2)

## The package contains the following parts:

1. Examination lamp equipped with halogen bulb
2. Screws for fixing (x 4)
3. User Manual.

## Uses:

This lamp is specially designed for indoor lighting, use it with our various fitting optional equipments.

## Recommendations:

Before installing the lamp or maintenance (for example: bulb replacement), always disconnect the plug from the socket.

Do not use the examination lamp in a humid atmosphere.

When bulb replacement please check the wattage: it must not exceed the authorized power.

Do not close ventilation openings on the head lamp.

Bulb replacement: wait until it has cooled sufficiently to avoid any risk of burns.

## Maintenance instructions and cleaning:

Mounting and connecting the lamp examination should be performed by qualified persons. For any maintenance, be sure that the lamp is unplugged.

All cleaning products (disinfectant, detergent, etc ....) can be used. For any rebel spot do not hesitate to use an abrasive cream agent.

## Installation instructions:

1. Pull out the plastic cover of the lamp base
2. Screw the lamp on the fixation (by the top)
3. Replace the plastic cover



If the external cable of this lamp is damaged, it must be replaced exclusively by the manufacturer, its service agent or similarly qualified person, in order to avoid any risk



If the external cable of this lamp is damaged, it must be replaced exclusively by the manufacturer, its service agent or similarly qualified person, in order to avoid any risk.



Review date : 05/2014